







Relative X-Ray Reflectometry for Characterization of Nanostructures

**A.G. Touriyanski, I.V. Pirshin,
N.L. Popov, Yu.A. Uspenski,
A.V. Vinogradov**



TWO-WAVE X-RAY REFLECTOMETER Making Reflectometry Versatile and Reliable (USA Patent 6041098/RF Patent 2124481)

COMPARISON OF EXISTING SCHEMES

The essence of the *Two-Wave Reflectometer (2-WXR)* scheme can be easily understood from comparison between the presently available *Philips X'Pert Reflectivity Diffractometer and Bruker Reflectometer* (Fig. 1a) and 2-WXR (Fig. 1b).

Both devices employ polychromatic radiation for the sample irradiation and separate characteristic lines from the reflected beam by PG monochromators. However, the *Philips/Bruker systems* use only one characteristic line, but the 2-WXR detects simultaneously two. This is achieved thanks to semitransparent HOPG monochromators and the new patented X-ray arrangement.

NEW MEASUREMENT OPPORTUNITIES

- 1. Double volume of information**
At equivalent geometry and sample irradiation conditions 2-WXR always gives twice as much information as conventional one (Fig 2).
- 2. Ratio reflectivity mode**
Intensity values and its ratio are measured simultaneously for two spectral lines at every angle step. It provides a complete suppression of the X-ray source instability and illumination variations with angle.
- 3. Refractometry mode.**
Direct determination of refractive index of surface layers by transmitting a probing beam through a sample edge. (Fig. 3).

NEW SCOPE OF APPLICATIONS

The above-mentioned 2-WXR features offer new opportunities in investigation methods and allow controlling samples that have been considered earlier as unacceptable for X-ray optical analysis.

- ◆ *Small samples*
- ◆ *Irregular form samples*
- ◆ *Locally processed samples*
- ◆ *Deeply buried or capsulated layers*
- ◆ *Curved samples*

CONCLUSION

Presented results clearly demonstrate that *the Two-Wave Reflectometer* and relevant measurement techniques offer a powerful tool for diagnostics of surfaces and thin film structures with *unmatched versatility, reliability, and information volume*.

Cont: e-mail: tour@mail1.lebedev.ru, t:(095)132-6268,f:(095)135-7880

Figures

Fig.1 a, b – X-ray schemes of the X'Pert Reflectivity Diffractometer from Philips Analytical X-Ray BV and Two-Wave Reflectometer respectively.

Fig. 2. Reflectivity curves from Ni-C multilayer on the Si substrate between 1st Bragg maximum and the critical angle of total external reflection for CuK_α and CuK_β lines (period d=5.16 nm, number of periods N=52).

Fig. 3. Refracted intensity vs. deviation angle (refractogram) for a C - Ni bilayer on Si substrate (C - 36 nm, Ni - 115 nm).

Fig. 4. Angle dependence of the reflectivity ratio R(CuK_α) / R(CuK_β) for a Si wafer fragment (5 x 7 mm²).

Fig. 5. Angle dependence of the reflectivity ratio R(CuK_α) / R(CuK_β) for an irregular form fragment of epitaxial Si_{1-x}Ge_x/Si structure (S≈10 mm²): dots - experiment, solid curve - calculations.

Fig. 6. Angle dependencies of the modified reflectivity ratio [R(β)-r(β)] / [R(α)-r(α)] for the ion-implanted diamond plate depicted in the upper right corner. R and r are respectively coefficients of reflection from an ion-implanted area and pure diamond; 1 - experiment, 2 - calculated.

Fig. 7. Refracted intensity vs. deviation angle from multilayer heterostructure Ge_xSi_{1-x}/Si grown by MBE for CuK_α and CuK_β lines.

Fig. 8a and 8b. Angle diagrams and corresponding PSD functions of the standard Si wafer for CuK_α and CuK_β lines characterizing surface roughness and radius of curvature).

Fig. 9a and 9b. Angle dependencies of reflectivity for CuK and CuK lines for the diamond-like C film (54 nm) on Si substrate and corresponding a reflectivity ratio R(CuK_α) / R(CuK_β) .

*P.N. Lebedev Physical Institute, 53,
Leninskii pr., Moscow, 119991, Russia,
tour@mail1.lebedev.ru*

